

UltrAi Viewer User Guide - Version: 0.3.0

Notes and Disclaimers

- This is a beta release; features and results may change.
- Results are for research use only and require clinical validation.
- This guide is mostly AI-generated and may have errors.

Supported OS

- Windows: [download](#)

macOS and Linux are currently not supported.

Compute Requirements

- CPU-only.

Supported files

- Image type: ultrasound images and videos.
- Preferred format: mp4, avi, TIFF, PNG, JPEG.

Quick Start

1. Install the App by double clicking the installer and following the instructions.
2. Launch the app just like any other Windows App.
3. Load a video: click on the **“Load video(s)”** button in the left panel.
4. Advance through the video using the playback options (slider, keyboard arrows, slider buttons)
5. Mark the start of a new clip by clicking the **“Start new clip”** button.
6. Click the **“Postprocess clips”** button to split the clips into their separate videos and classes. Results will be saved in the selected output folder.

Other features

7. Segment the target object frame by frame using the mouse. Note: the masks will be propagated to the next frame.
8. Edit the annotations if needed using the Edit tools: Freehand Line, Segmented Line, Paint Brush, Eraser.
9. Save the mask: click the **“Save masks”** button. Note, **results** are **saved automatically** to the output folder as you advance through the video.

Editing Tools

- Freehand Line: draw an outline ROI.
- Segmented Line: click to add points, right-click or double-click to finish.
- Paint Brush: add mask.

- Eraser: remove mask.
- Tool Radius: controls edit thickness.
- Fill ROI: unchecked shows contour; checked fills the ROI into the mask.
- Ctrl+Z to undo, Ctrl+Y to redo.